

SECTION 1: Identification of the substance/mixture and of the company/undertaking**Identification of the substance or mixture**

Product name Master Mix

Company/undertaking identification

Applied DNA Sciences
50 Health Sciences Drive
Stony Brook, NY 11790
Phone: +1-631-240-8800
Email: safety@adnas.com
<http://www.adnas.com/>

**24 hour Emergency Response for Hazardous Materials
[or Dangerous Goods] Incident. Spill, Leak, Fire,
Exposure, or Accident.**

631-240-8800 (Monday through Friday,
from 8:30 A.M. to 5:00 P.M. EST)

SECTION 2: Hazards identification**GHS Classification****Signal Word**

None

Hazard pictograms

None

Health hazards

Not classified

Physical hazards

Not classified

Environmental hazards

Not classified

Hazard Statements

Not Applicable

Precautionary Statements**Prevention**

Not Applicable

Response

Not Applicable

Storage

Not Applicable

Disposal

Not Applicable

Other hazards

Not Applicable

HMIS

Health	0
Flammability	0
Reactivity	0

SECTION 3: Composition/information on ingredients

Component	CAS No.	Common name	EINECS-No	Weight-%
Glycerol 56-81-5 (10-30)	56-81-5	-	200-289-5	10-30

We recommend handling all chemicals with caution.

SECTION 4: First aid measures**Description of first aid measures**

Skin contact	Rinse skin with water. Immediate medical attention is not required.
Eye contact	Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
Ingestion	Not expected to present a significant ingestion hazard under anticipated conditions of normal use. If you feel unwell, seek medical advice.
Inhalation	Not expected to be an inhalation hazard under anticipated conditions of normal use of this material. Consult a physician if necessary.
Notes to Physician	Treat symptomatically.

Most important symptoms and effects, both acute and delayed

Not Applicable

Indication of any immediate medical attention and special treatment needed

None.

SECTION 5: Firefighting measures**Extinguishing media**

Suitable extinguishing media
Unsuitable extinguishing media

Water spray. Carbon dioxide (CO₂). Foam. Dry chemical.
No information available.

Special hazards arising from the substance or mixture

Not known

Protective equipment and precautions for firefighters

Standard procedure for chemical fires.

SECTION 6: Accidental release measures**Personal precautions, protective equipment and emergency procedures**

Ensure adequate ventilation
Always wear recommended Personal Protective Equipment
Use personal protection equipment
See section 8 for more information

Environmental precautions

No special environmental precautions required.

Methods and material for containment and cleaning up

Soak up with inert absorbent material.

Reference to other sections

See section 8 for more information.

SECTION 7: Handling and storage

Precautions for safe handling

Use personal protective equipment as required. No special handling advices are necessary.

Conditions for safe storage, including any incompatibilities

Keep in a dry, cool and well-ventilated place. Keep in properly labeled containers.

Storage Conditions

Store at -25 C to -15 C.

SECTION 8: Exposure controls/personal protection

Control parameters

Chemical Name	OSHA PEL	OSHA PEL (Ceiling)	ACGIH OEL (TWA)	ACGIH OEL (STEL)
Glycerol	15 mg/m ³ 5 mg/m ³	None	None	None

Chemical Name	Brazil - OEL - TWAs (LTs)	Brazil - OEL - Ceilings	Brazil - OEL - Skin Designations
Glycerol	None	None	None

Engineering Measures Ensure adequate ventilation, especially in confined areas.

Exposure controls
Personal Protective Equipment

Respiratory protection In case of insufficient ventilation wear respirators and components tested and approved under appropriate government standards.

Hand protection Wear suitable gloves. Glove material: Compatible chemical-resistant gloves.

Eye protection Tight sealing safety goggles.

Skin and Body Protection Wear suitable protective clothing.

Hygiene Measures Handle in accordance with good industrial hygiene and safety practice

Environmental exposure controls

No special environmental precautions required.

SECTION 9: Physical and chemical properties**Information on basic physical and chemical properties**

Appearance	liquid	
Odor	No data available	
Odor Threshold	No data	
Molecular Weight	No data	
Melting point / melting range	°C No data	°F No data
Boiling point / boiling range	°C No data	°F No data
Flash point	°C No data	°F No data
Autoignition Temperature	°C No data	°F No data
Decomposition temperature	°C No data	°F No data
Evaporation rate	No data	
Flammability (solid, gas)	No data available	
Upper explosion limit	No data	
Lower explosion limit	No data	
Vapor Pressure	No data	
Vapor density	No data	
Relative density	No data	
Specific gravity	No data	
Solubility	No data available	
Partition coefficient: n-octanol/water	No data available	
Viscosity	No data	
Explosive properties	No data	
Oxidizing properties	No data	

Other information

No data available.

SECTION 10: Stability and reactivity

Reactivity	None known.
Chemical stability	Stable under normal conditions.
Possibility of hazardous reactions	Hazardous reaction has not been reported.
Conditions to avoid	No information available.
Incompatible materials	No dangerous reaction known under conditions of normal use.
Hazardous decomposition products	No data available.

SECTION 11: Toxicological information
Information on toxicological effects

Chemical Name	Oral LD50	Dermal LD50	Inhalation LC50
Glycerol	= 12600 mg/kg (Rat)	No data available	>570mg/m ³ (Rat)

Principal Routes of Exposure

Acute toxicity Data are conclusive but insufficient for classification.

Skin corrosion/irritation Data are conclusive but insufficient for classification

Serious eye damage/irritation Data are conclusive but insufficient for classification

Respiratory or skin sensitization Data are conclusive but insufficient for classification

Specific target organ toxicity (STOT) – single exposure Data are conclusive but insufficient for classification

Specific target organ toxicity (STOT) – repeated exposure Data are conclusive but insufficient for classification

Carcinogenicity Data are conclusive but insufficient for classification

Germ cell mutagenicity Data are conclusive but insufficient for classification

Reproductive toxicity Data are conclusive but insufficient for classification

Aspiration hazard Data are conclusive but insufficient for classification

SECTION 12: Ecological information
Ecotoxicity

Contains no substances known to be hazardous to the environment or not degradable in waste water treatment plants.

Chemical Name	Toxicity to algae	Toxicity to daphnia and other aquatic invertebrates	Toxicity to fish	Microtox Data	log Pow
Glycerol	No data available	Daphnia magna EC50>500 mg/L (24 h)	No data available	No data available	logPow-1.76

Mobility in soil No information available.

Persistence and degradability No information available.

Bioaccumulative potential No information available.

Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be a PBT or a vPvB.

Other adverse effects No information available.

SECTION 13: Disposal considerations
Waste treatment methods

The generation of waste should be avoided or minimized wherever possible. Empty containers or liners may retain some product residues. This material and its container must be disposed of in accordance with approved disposal technique. Disposal of this product, its solutions or of any by-products, shall comply with the requirements of all applicable local, regional or national/federal regulations

SECTION 14: Transport information
IATA / ADR / DOT-US / IMDG

Not regulated in the meaning of transport regulations

UN number	Not Applicable
UN proper shipping name	Not Applicable
Transport hazard class(es)	Not Applicable
Packing group	Not Applicable

Environmental hazards

Not Applicable

Special precautions for user

Not Applicable

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not Applicable.

SECTION 15: Regulatory information

Component	US TSCA
Glycerol 56-81-5 (10-30)	Listed

US Federal Regulations
SARA 313

This product is not regulated by SARA.

Clean Air Act, Section 112 Hazardous Air Pollutants (HAPs) (see 40 CFR 61)

This product does not contain HAPs

US State Regulations

California Proposition 65

This product does not contain any Proposition 65 chemicals.

WHMIS Hazard Class

Non-controlled

This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations (CPR) and the SDS contains all the information required by the CPR.

National Regulations - Brazil

<u>Chemical Name</u>	<u>CAS No.</u>	Brazil - National Agency for Sanitary Surveillance (ANVISA)	Brazil - National List of Carcinogen Agents to Humans (LINACH)
Glycerol	56-81-5	Not Listed	Not Listed

SECTION 16: Other information

Reason for revision New

Revision date 14-May-2020

References

- ECHA: <http://echa.europa.eu/>
- TOXNET: <http://toxnet.nlm.nih.gov/>
- eChemPortal: <http://www.echemportal.org/>
- LOLI database: <https://www.chemadvisor.com/loli-database>

"The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS SDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE"

End of Safety Data Sheet